

PRODUCT INFORMATION

Uniprot ID	P28907
Common Name	HuMax-CD38
Conjugate	PE-conjugated
Synonyms	CD38;T10;ADPRC1
Applications	Flow Cyt
Recommended Dilutions	1µl/test
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	IgG1
Reactivity	Human
Target	CD38
Description	PE-conjugated Anti-CD38 (daratumumab biosimilar) mAb
Delivery	In Stock
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.






Figure 1. Flow cytometry analysis with 1 μ l/test PE-conjugated Anti-CD38 (daratumumab biosimilar) on H929 cells(Blue histogram) or CHO-S cells (Red histogram).

DIMABIO CONFIDENTIAL

