

PRODUCT INFORMATION

| | |
|---|---|
| Uniprot ID | P28907 |
| Common Name | MOR_202;MOR202;MOR-202;MOR-03087;MOR 03087;MOR03087;MOR 3087;MOR-3087;MOR3087;TJ-202;TJ202;TJ 202 |
| Synonyms | CD38;ADPRC1;cADPR1;ADPRC 1;T10;cADPr hydrolase 1 |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid□PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Whole mAb |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | CD38 |
| Description | PE-conjugated Anti-CD38(felzartamab biosimilar) mAb |
| Delivery | 6-8 weeks |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | PE-conjugated |

