

PRODUCT INFORMATION

Uniprot ID	P49961
Common Name	ES 002;ES002;ES-002;ES 002023;ES-002023;ES 002023
Synonyms	CD39;SPG64;ATP-DPH;ATPDase;NTPDase-1;ENTPD1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CD39
Description	PE-conjugated Anti-CD39(eurestobart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

