

PRODUCT INFORMATION

| | |
|---|---|
| Uniprot ID | P25942 |
| Common Name | CFZ-533, NVP-CFZ533, OM11-62MF |
| Synonyms | CD40;Bp50;CDW40;MGC9013;TNFRSF5;p50 |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | N/A |
| Host Species | Homo sapiens |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | CD40 |
| Description | PE-conjugated Anti-CD40 (iscalimab biosimilar) mAb |
| Delivery | 3-4 weeks |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | PE-conjugated |

DIMABIO CONFIDENTIAL

