

PRODUCT INFORMATION

Uniprot ID	P29965
Common Name	INX 021;INX021;INX-021;SAR-441344;SAR441344;SAR 441344
Synonyms	IGM;IMD3;TRAP;gp39;CD154;CD40L;HIGM1;T-BAM;TNFSF5;hCD40L
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CD40 Ligand
Description	PE-conjugated Anti-CD40 Ligand(frexalimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

