

PRODUCT INFORMATION

Uniprot ID	P25942
Common Name	LVGN-7409;LVGN7409
Synonyms	CD40;p50;Bp50;CDW40;TNFRSF5
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CD40
Description	PE-conjugated Anti-CD40(dalnicastobart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

