

**PRODUCT INFORMATION**

<b>Uniprot ID</b>	P25942
<b>Common Name</b>	APX 005;APX005;APX 005;APX5;APX-5;APX 5;EPI-0050;EPI0050;EPI 0050;EPI 50;EPI-50;EPI50
<b>Synonyms</b>	CD40;p50;Bp50;CDW40;TNFRSF5
<b>Applications</b>	N/A
<b>Recommended Dilutions</b>	N/A
<b>Formulation &amp; Reconstitution</b>	N/A
<b>Host Species</b>	Whole mAb
<b>IgG type</b>	IgG1
<b>Reactivity</b>	Human
<b>Target</b>	CD40
<b>Description</b>	PE-conjugated Anti-CD40(sotigalimab biosimilar) mAb
<b>Delivery</b>	6-8 weeks
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated

DIMABIO CONFIDENTIAL

