

PRODUCT INFORMATION

Uniprot ID	Q08722
Common Name	TJ-011133, TJC-4
Synonyms	CD47;MER6;IAP
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	N/A
Host Species	Humanized
IgG type	IgG4
Reactivity	Human
Target	CD47
Description	PE-conjugated Anti-CD47(Iemzoparlimab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

