

**PRODUCT INFORMATION**

<b>Clone ID</b>	DM161
<b>Target</b>	CD5
<b>Synonyms</b>	CD5;LEU1
<b>Host Species</b>	Rabbit
<b>Description</b>	PE-conjugated Anti-CD5 antibody(DM161); Rabbit mAb
<b>Delivery</b>	3-4 weeks
<b>Uniprot ID</b>	P06127
<b>IgG type</b>	Rabbit IgG
<b>Clonality</b>	Monoclonal
<b>Reactivity</b>	Human
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Purification</b>	Purified from cell culture supernatant by affinity chromatography
<b>Formulation &amp; Reconstitution</b>	Liquid□PBS with 0.05% Proclin300, 1% BSA
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	<p>This gene encodes a member of the scavenger receptor cysteine-rich (SRCR) superfamily. Members of this family are secreted or membrane-anchored proteins mainly found in cells associated with the immune system. This protein is a type-I transmembrane glycoprotein found on the surface of thymocytes; T lymphocytes and a subset of B lymphocytes. The encoded protein contains three SRCR domains and may act as a receptor to regulate T-cell proliferation. Alternative splicing results in multiple transcript variants encoding different isoforms.</p>
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated

