

PRODUCT INFORMATION

Uniprot ID	P32970
Common Name	ARGX-110;ARGX110;JNJ 74494550;JNJ74494550;JNJ-74494550;JNJ-4550;JNJ4550
Synonyms	CD27L;LPFS3;CD27-L;CD27LG;TNFSF7;TNLG8A
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CD70
Description	PE-conjugated Anti-CD70(cusatuzumab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

