

PRODUCT INFORMATION

Clone ID	DM204
Target	CD73
Synonyms	CD73;NT5E;5'-Nucleotidase;5'-NT;NT5;NTE
Host Species	Rabbit
Description	PE-conjugated Anti-CD73 antibody(DM204); Rabbit mAb
Delivery	3-4 weeks
Uniprot ID	P21589
IgG type	Rabbit IgG
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	The protein encoded by this gene is a plasma membrane protein that catalyzes the conversion of extracellular nucleotides to membrane-permeable nucleosides. The encoded protein is used as a determinant of lymphocyte differentiation. Defects in this gene can lead to the calcification of joints and arteries. Two transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq; Mar 2011]
Usage	Research use only
Conjugate	PE-conjugated

