

PRODUCT INFORMATION

Uniprot ID	P21589
Common Name	COR-004;COR004;COR 004;COR-4;COR4;COR 4;CPI-006;CPI006;CPI 006;CPI-6;CPI6;CPI 6;CPX-006;CPX006;CPX 006;CPX-6;CPX6;CPX 6
Synonyms	NT;eN;NT5;NTE;eNT;CD73;E5NT;CALJA
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CD73
Description	PE-conjugated Anti-CD73(mupadolimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

