

**PRODUCT INFORMATION**

<b>Uniprot ID</b>	P21589
<b>Common Name</b>	MEDI-9447, MEDI9447
<b>Synonyms</b>	NT5E, 5'-Nucleotidase, 5'-NT, NT5, NTE
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	Liquid □ PBS with 0.05% Proclin300, 1% BSA
<b>Host Species</b>	Homo sapiens
<b>IgG type</b>	IgG1(K97R,L117F,L118E,P214S)
<b>Reactivity</b>	Human
<b>Target</b>	CD73
<b>Description</b>	PE-conjugated Anti-CD73(oleclumab biosimilar) mAb
<b>Delivery</b>	3-4 weeks
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated

