

PRODUCT INFORMATION

Uniprot ID	P04233
Common Name	CD74-DOX (ADC),hLL1,hLL1-DOX (ADC),MEDI-115
Synonyms	DHLAG
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1(K97R)
Reactivity	Human
Target	CD74
Description	PE-conjugated Anti-CD74(milatumzumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

