

PRODUCT INFORMATION

Uniprot ID	P40199
Common Name	BAY-1834942
Synonyms	CD66c;CEAL;NCA
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG2
Reactivity	Human
Target	CEACAM6
Description	PE-conjugated Anti-CEACAM6(tinurilimab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

