

PRODUCT INFORMATION

Clone ID	DMC442
Target	CHODL
Synonyms	C21orf68; MT75; PRED12
Host Species	Rabbit
Description	PE-conjugated Anti-CHODL antibody(DMC442); IgG1 Chimeric mAb
Delivery	3-4 weeks
Uniprot ID	Q9H9P2
IgG type	Rabbit/Human Fc chimeric IgG1
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	This gene encodes a type I membrane protein with a carbohydrate recognition domain characteristic of C-type lectins in its extracellular portion. In other proteins; this domain is involved in endocytosis of glycoproteins and exogenous sugar-bearing pathogens. This protein localizes predominantly to the perinuclear region. Several transcript variants encoding a few different isoforms have been found for this gene.
Usage	Research use only
Conjugate	PE-conjugated

