

PRODUCT INFORMATION

Uniprot ID	P56747
Common Name	ASP1650, IMAB027
Synonyms	Claudin 6;Claudin-6;Skullin
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	Human IgG1
Reactivity	Human
Target	CLDN6
Description	PE-conjugated Anti-CLDN6 (IMAB027) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMA BIO CONFIDENTIAL

