

PRODUCT INFORMATION

Uniprot ID	O95727
Common Name	HuC110-34
Synonyms	CD355
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1 (L115E)
Reactivity	Human
Target	CRTAM
Description	PE-conjugated Anti-CRTAM(HuC110-34 biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

