

PRODUCT INFORMATION

Uniprot ID	P16410
Common Name	PSB105;PSB-105;PSB-205 (cocktail with Tuvonralimab);PSB205 (cocktail with Tuvonralimab);QL1706 (cocktail with Tuvonralimab)
Synonyms	CD;GSE;GRD4;ALP55;CD152;CTLA-4;IDDM12;CELIAC3
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CTLA-4
Description	PE-conjugated Anti-CTLA-4(tuvonralimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

