

PRODUCT INFORMATION

Uniprot ID	P00533
Common Name	IMC-11F8;IMC 11F8;IMC11F8;LY3012211;LY3012211;LY-3012211;Portraza;Portraza
Synonyms	ERBB;ERRP;HER1;mENA;ERBB1;NNCIS;PIG61;NISBD2
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	N/A
IgG type	IgG1
Reactivity	Human
Target	EGFR
Description	PE-conjugated Anti-EGFR(necitumumab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

