

PRODUCT INFORMATION

Clone ID	DM154
Target	FAP
Synonyms	FAP;FAPalpha;SIMP;Seprase;APCE
Host Species	Rabbit
Description	PE-conjugated Anti-FAP antibody(DM154); Rabbit mAb
Delivery	3-4 weeks
Uniprot ID	Q12884
IgG type	Rabbit IgG
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	<p>The protein encoded by this gene is a homodimeric integral membrane gelatinase belonging to the serine protease family. It is selectively expressed in reactive stromal fibroblasts of epithelial cancers; granulation tissue of healing wounds; and malignant cells of bone and soft tissue sarcomas. This protein is thought to be involved in the control of fibroblast growth or epithelial-mesenchymal interactions during development; tissue repair; and epithelial carcinogenesis. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.</p>
Usage	Research use only
Conjugate	PE-conjugated

