

PRODUCT INFORMATION

| | |
|---|---|
| Clone ID | DMC220 |
| Target | FCRL5 |
| Synonyms | BXMAS1; CD307; CD307e; FCRH5; IRTA2; PRO820 |
| Host Species | Rabbit |
| Description | PE-conjugated Anti-FCRL5 antibody(DMC220); IgG1 Chimeric mAb |
| Delivery | Under Development |
| Uniprot ID | Q96RD9 |
| IgG type | Rabbit/Human Fc chimeric IgG1 |
| Clonality | Monoclonal |
| Reactivity | Human |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Purification | Purified from cell culture supernatant by affinity chromatography |
| Formulation & Reconstitution | Liquid□PBS with 0.05% Proclin300, 1% BSA |
| Storage & Shipping | Store at 2°C-8°C for 6 months |

Background This gene encodes a member of the immunoglobulin receptor superfamily and the Fc-receptor like family. This gene and several other Fc receptor-like gene members are clustered on the long arm of chromosome 1. The encoded protein is a single-pass type I membrane protein and contains 8 immunoglobulin-like C2-type domains. This gene is implicated in B cell development and lymphomagenesis. Alternatively spliced transcript variants encoding different isoforms have been identified.

Usage Research use only

Conjugate PE-conjugated

DIMA Disclaimer All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

