

PRODUCT INFORMATION

Uniprot ID	Q96RD9
Common Name	BFCR-4350-A, BFCR-4350A, RG-6160, RO-7187797
Synonyms	BXMAS1;CD307;CD307e;FCRH5;IRTA2;PRO820
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	FCRL5
Description	PE-conjugated Anti-FCRL5(cevostamab biosimilar without anti-CD3 portion) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

