

PRODUCT INFORMATION

Uniprot ID	P02751
Common Name	PYX-201;PYX201;PYX 201
Synonyms	FN;CIG;FNZ;MSF;ED-B;FINC;GFND;LETS;GFND2;SMDCF
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	FN1
Description	PE-conjugated Anti-FN1(micvotabart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

