

PRODUCT INFORMATION

Uniprot ID	Q99988
Common Name	AV-380;AV 380;AV380
Synonyms	HG;PDF;MIC1;PLAB;MIC-1;NAG-1;PTGFB;GDF-15
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	GDF15
Description	PE-conjugated Anti-GDF15(rilogrotug biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

