

PRODUCT INFORMATION

Clone ID	DMC478
Target	GPR87
Synonyms	FKSG78; GPR95; KPG_002
Host Species	Rabbit
Description	PE-conjugated Anti-GPR87 antibody(DMC478); IgG1 Chimeric mAb
Delivery	3-4 weeks
Uniprot ID	Q9BY21
IgG type	Rabbit/Human Fc chimeric IgG1
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	This gene encodes a G protein-coupled receptor and is located in a cluster of G protein-coupled receptor genes on chromosome 3. The encoded protein has been shown to be overexpressed in lung squamous cell carcinoma (PMID:18057535) and regulated by p53 (PMID:19602589). [provided by RefSeq; Nov 2011]
Usage	Research use only
Conjugate	PE-conjugated

