

PRODUCT INFORMATION

Uniprot ID	Q9NZD1
Common Name	LBL-034
Synonyms	GPRC5D
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1(L117A, L118A, D148A, P212A, Y232C, E239D, M241L, T249S, L251A, Y290V)
Reactivity	Human
Target	GPRC5D
Description	PE-conjugated Anti-GPRC5D(LBL-034 biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

