

**PRODUCT INFORMATION**

<b>Uniprot ID</b>	P08069
<b>Common Name</b>	RO4858696-000,teprotumumab-trbw
<b>Synonyms</b>	CD221, IGFIR, IGFR, JTK13
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	Liquid□PBS with 0.05% Proclin300, 1% BSA
<b>Host Species</b>	Homo sapiens
<b>IgG type</b>	IgG1(E239D, M241L)
<b>Reactivity</b>	Human
<b>Target</b>	IGF-1R
<b>Description</b>	PE-conjugated Anti-IGF-1R(teprotumumab biosimilar) mAb
<b>Delivery</b>	3-4 weeks
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated

