

PRODUCT INFORMATION

Uniprot ID	P35225
Common Name	ABT-308;ABT308;CC 93538;CC93538;CC-93538;RPC 4046;RPC4046;RPC-4046
Synonyms	P600;IL-13;NC30
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	IL13
Description	PE-conjugated Anti-IL13(cendakimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

