

PRODUCT INFORMATION

Uniprot ID	P35225
Common Name	DRM 06;DRM-06;DRM06;DRM 6;DRM-6;DRM6;LY 3650150;LY-3650150;LY3650150;MILR1444A;MILR-1444A;MILR 1444;ARG 3637;ARG-3637;ARG3637;RO 5490255;RO-5490255;RO5490255;TNX 650;TNX-650;TXN650
Synonyms	P600;IL-13;NC30
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	IL13
Description	PE-conjugated Anti-IL13(lebrikizumab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

