

PRODUCT INFORMATION

Uniprot ID	Q16552
Common Name	AK-111;AK111;AK 111
Synonyms	IL17;CTLA8;IL-17;ILA17;CTLA-8;IL-17A
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	IL17A
Description	PE-conjugated Anti-IL17A(gumokimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

