

PRODUCT INFORMATION

Uniprot ID	Q96F46
Common Name	AMG827, KHK4827
Synonyms	CD217;CDw217;IL-17RA;IL17R;CANDF5;hIL-17R
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	IgG2
Reactivity	Human
Target	IL17RA
Description	PE-conjugated Anti-IL17RA (brodalumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

