

PRODUCT INFORMATION

Uniprot ID	Q9HB29
Common Name	BI-655130, BI655130, BI 655130, spesolimab-sbzo
Synonyms	ILRL2; IL-36R; IL1RRP2; IL-36R
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid □ PBS with 0.05% Proclin300, 1% BSA
Host Species	Chimeric Humanized
IgG type	IgG1(K97R,L117A,L118A)
Reactivity	Human
Target	IL1RL2
Description	PE-conjugated Anti-IL1RL2(spesolimab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

