

PRODUCT INFORMATION

Uniprot ID	Q8N6P7
Common Name	LEO 138559;LEO-138559;LEO138559;ARGX112;ARGX-112;LP0145;LP-0145
Synonyms	IL22R;CRF2-9;IL22R1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	IL22RA1
Description	PE-conjugated Anti-IL22RA1(temtokibart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

