

PRODUCT INFORMATION

Uniprot ID	Q9NPF7
Common Name	CNTO-1959
Synonyms	SGRF
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	N/A
Host Species	Humanized
IgG type	IgG1(E239D, M241L)
Reactivity	Human
Target	IL23A
Description	PE-conjugated Anti-IL23A(guselkumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

