

PRODUCT INFORMATION

Uniprot ID	Q9NPF7
Common Name	IBI-112;IBI112;IBI 112
Synonyms	P19;SGRF;IL-23;IL-23A;IL23P19
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	IL23A
Description	PE-conjugated Anti-IL23A(picankibart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

