

PRODUCT INFORMATION

Uniprot ID	O95760
Common Name	REGN-3500;REGN3500;REGN 3500;SAR-440340;SAR440340;SAR 440340
Synonyms	DVS27;IL1F11;NF-HEV;NFEHEV;C9orf26
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	IL33
Description	PE-conjugated Anti-IL33(itepekimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

