

PRODUCT INFORMATION

Uniprot ID	P05113
Common Name	AQ82742999;AQ-82742999;GSK 294;GSK-294;GSK294;GSK-3511294;GSK3511294
Synonyms	EDF;TRF;IL-5;Interleukin-5
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	IL5
Description	PE-conjugated Anti-IL5(depemokimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

