

PRODUCT INFORMATION

Uniprot ID	Q01344
Common Name	BIW-8405,KHK4563,MEDI-563
Synonyms	CDw125; CD125; IL5R
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1(E239D,M241L)
Reactivity	Human
Target	IL5RA
Description	PE-conjugated Anti-IL5RA(benralizumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

