

PRODUCT INFORMATION

Uniprot ID	P16871
Common Name	ZB-168;ZB168;ZB 168
Synonyms	ILRA;CD127;IL7RA;CDW127;IMD104;sIL-7R;Inc-IL7R;IL7Ralpha;IL-7Ralpha;IL-7R-alpha
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	IL7RA
Description	PE-conjugated Anti-IL7RA(crebankitug biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

