

PRODUCT INFORMATION

Uniprot ID	Q71H61
Common Name	BAY-1905254
Synonyms	C1orf32; Angulin-3
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	IgG2(S257A)
Reactivity	Human
Target	ILDR2
Description	PE-conjugated Anti-ILDR2(bapotulimab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

