

PRODUCT INFORMATION

Uniprot ID	P18627
Common Name	LAG-3,CD223
Synonyms	BMS-986016;ONO-4482
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	IgG4
Reactivity	Human
Target	LAG3
Description	PE-conjugated Anti-LAG3 (relatlimab biosimilar)mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

