

PRODUCT INFORMATION

Uniprot ID	P18627
Common Name	BSI-045B;BSI045B;BSI 045B;TQC-2731;TQC2731;TQC 2731
Synonyms	LAG-3;CD223
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	LAG3
Description	PE-conjugated Anti-LAG3(brivestobart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

