

**PRODUCT INFORMATION**

<b>Uniprot ID</b>	P18627
<b>Common Name</b>	BSI-045B;BSI045B;BSI045B;TQC-2731;TQC2731;TQC 2731
<b>Synonyms</b>	LAG-3;CD223
<b>Applications</b>	N/A
<b>Recommended Dilutions</b>	N/A
<b>Formulation &amp; Reconstitution</b>	N/A
<b>Host Species</b>	Whole mAb
<b>IgG type</b>	IgG4
<b>Reactivity</b>	Human
<b>Target</b>	LAG3
<b>Description</b>	PE-conjugated Anti-LAG3(brivestobart biosimilar) mAb
<b>Delivery</b>	6-8 weeks
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only

DIMABIO CONFIDENTIAL

