

PRODUCT INFORMATION

Uniprot ID	P18627
Common Name	IMP-701;IMP 701;IMP701;LAG-525;LAG525;LAG 525;NVP-LAG525
Synonyms	LAG-3;CD223
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	LAG3
Description	PE-conjugated Anti-LAG3(ieramilimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

