

PRODUCT INFORMATION

Clone ID	DMC277
Target	LGALS1
Synonyms	GAL1; GBP
Host Species	Rabbit
Description	PE-conjugated Anti-LGALS1 antibody(DMC277); IgG1 Chimeric mAb
Delivery	3-4 weeks
Uniprot ID	P09382
IgG type	Rabbit/Human Fc chimeric IgG1
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. This gene product may act as an autocrine negative growth factor that regulates cell proliferation.
Usage	Research use only
Conjugate	PE-conjugated

