

PRODUCT INFORMATION

Uniprot ID	Q8NHJ6
Common Name	MK-0482-Hu52B8
Synonyms	ILT3; LIR5; CD85k
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid-PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG4
Reactivity	Human
Target	LILRB4
Description	PE-conjugated Anti-LILRB4(MK-0482 biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

