

PRODUCT INFORMATION

Uniprot ID	Q16553
Common Name	ANTI-LY6E□DLYE-5953A□DLYE5953A□RG 7841
Synonyms	RIGE; SCA2; RIG-E; SCA-2; TSA-1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	LY6E
Description	PE-conjugated Anti-LY6E(RG7841 biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

