

PRODUCT INFORMATION

Uniprot ID	Q04912
Common Name	IMC RON-8, IMC-RON8, RON8
Synonyms	CD136;CDw136;NPCA3;PTK8;RON
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	IgG1
Reactivity	Human
Target	MST1R
Description	PE-conjugated Anti-MST1R(narnatumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

