

PRODUCT INFORMATION

Uniprot ID	P26715
Common Name	IPH-2201, IPH2201, NN-8765, NNC-0141-0000-0100, NNC-141-01000, anti-NKG2A, humZ270
Synonyms	CD159A;NKG2;NKG2A
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG4
Reactivity	Human
Target	NKG2A
Description	PE-conjugated Anti-NKG2A(monalizumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

